

Ideal® Veterinary Apparel

Designed for veterinary use when performing cesarean sections and other large animal procedures, Ideal veterinary apparel is made from a urethane coated nylon oxford material. The garments are navy blue in color, highly water resistant, breathable and machine washable. Custom sizes of all garments may be special ordered. Please specify size when ordering, by adding S, M, L or XL after the item number (for example, VA301-L).

Apparel Measurements

Tops

Size	Chest	Arms
SMALL	46"	18"
MEDIUM	49"	18"
LARGE	54"	18"
X-LARGE	55"	20"

Bottoms

Size	Waist	Inseam
MEDIUM	42"	31"
LARGE	44"	32"
X-LARGE	48"	34"

Cleaning Instructions

- Shake out after use
- Remove any excess dirt and grime
- Spray with hose and pre-soak
- Wash - normal cycle in COLD water
- Shake out - DO NOT DRY IN DRYER - line dry only

OB Jacket

Item No.	Description
VA301	OB Jacket (Please specify size when ordering: L or XL)

OB Suit

- Pants have an elastic waistband

Item No.	Description
VA101	Top only - specify size: M, L or XL
VA102	Bottoms only - specify size: M, L or XL

Vet Apron

- Extra long ties

Item No.	Description
VA400	Vet Apron - One size

OB Gown

- Velcro fasteners in back

Item No.	Description
VA500	OB Gown (Please specify size when ordering: L or XL)

Shoulder Protector

- Buckle with adjustable strap
- May be worn on either arm

Item No.	Description
VA900	Shoulder Protector (Please specify size when ordering: S or L)

Vet Overall

- Velcro fasteners on straps

Size	Waist	Inseam
LARGE	41"	33"
X-LARGE	47"	33½"

Item No.	Description
VA600	Vet Overall (Please specify size when ordering: L or XL)



Disposable Apparel

PolyGown™

Manufactured of 1.3 mil lightweight embossed polyethylene film, the Ideal PolyGown is designed for veterinary and livestock applications where minimal protection is needed. PolyGown is packaged in a convenient dispenser box for easy removal and is available in three sizes: M, L and XL.

Item No.	Description
DG100-M	PolyGown, size M - box of 20
DG100-L	PolyGown, size L - box of 15
DG100-XL	PolyGown, size XL - box of 10

Veterinary Gown

The Ideal Veterinary Gown is manufactured of durable 1.2 mil polyethylene film and is designed for those veterinary procedures where disposable full-body protection is needed. The Veterinary Gown is 61" in length, has a slit neck for ease in donning and comfort pleats in the back to facilitate ease of movement. The Veterinary Gown is packaged in a heat-sealed bag containing 25 gowns. One size.

Item No.	Description
VG100	Veterinary Gown, one size - pack of 25

Disposable Overalls

For use in contaminated areas, Ideal's non-woven disposable coveralls are available in white with a zipper front.

Item No.	Description
TA713	Medium
TA714	Large
TA715	X-Large
TA715XXL	XX-Large

PolyApron™

Designed for durable protection of clothing when treating or grooming large and small animals, the disposable Ideal PolyApron is made from blue heavyweight 4.5 mil polyethylene film. Available in one size (28" x 46"), the apron has extra-long ties for ease of tying in front. PolyApron is packaged and sealed in a polyethylene sleeve containing 10 aprons.

Item No.	Description
A450	PolyApron, one size - pack of 10

Blue Dairy Apron

Ideal's Blue Dairy Apron is available in blue with adjustable neck and back ties. One size fits all.

Item No.	Description
TA1858	Blue Dairy Apron

