

Spectrasol[®] Disinfectant



CLEANER, DISINFECTANT, DETERGENT, DEODORIZER, VIRUCIDE, FUNGICIDE (against pathogenic fungi), MILDEWSTAT (on hard, non-porous, inanimate surfaces)

A Phosphate-Free Quaternary Ammonium Formulation Designed To Provide Effective Cleaning, Deodorizing And Disinfection

Product

This tech sheet is intended for use only as a guide in providing general information regarding the directions, warnings and cautions associated with the use of this product. As with any disinfectant, always read and understand the entire label instructions on the package and MSDS before using.

Apply product to soiled area with a mop, cloth, sponge, or hand pump trigger sprayer. Thoroughly wet soiled surface and scrub as necessary. Rinse surface with potable water.

For heavily soiled areas, a pre-cleaning step is required. Prepare a fresh solution for each use. Apply solution with a mop, cloth, sponge, or hand pump trigger sprayer so as to wet all surfaces thoroughly. Allow to remain wet for 10 minutes, then remove excess liquid.

Active Ingredients

Octyl decyl dimethyl ammonium chloride	3.255%
Dioctyl dimethyl ammonium chloride	1.628%
Didecyl dimethyl ammonium chloride	1.628%
Alkyl (C ₁₄ -50%, C ₁₂ -40%, C ₁₆ -10%) dimethyl benzyl ammonium chloride	4.339%
Other ingredients:	89.150%
Total.....	100.000%

Physical Properties

State	Liquid
Color	Clear
Odor	Mild Mint
pH, 25°C	6.6
Density	8.41 lbs/gal; 1 kg/L
Boiling Point	200°F/93.3°C
Volatility	80%
Solubility	100%

EPA Reg. No. 1839-166-134

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Benefits

- Economical to use; one ounce of concentrate diluted in one gallon of water provides an effective solution for disinfecting.
- Cleans as it disinfects, and deodorizes areas that are prone to odors caused by microorganisms.
- Effective as a fungicide, virucide and mildewstat.
- Non-corrosive and safe to use on cages and tables; will not dull high-gloss floor finishes, even with repeated use.

Packaging

Packaging	Shipping Weight	Product Code
4 x 1 gal (3.79 L)/case	34 lb (15.4 kg)	190421

Usage Rates

Task	U.S.	Metric
General Use	1 oz/gal	8 mL/L

