



Acid Tray Wash EVO

ACIDIFIED, LIQUID TRAY WASH

Product

Acid Tray Wash EVO is an acidified, phosphate-free, non-foaming cleaner for washing trays, plastic chick boxes, egg cases, and filler flats in automatic tray washers. It is also an excellent acidic rinse and equipment descaler.

Usage Rates

Task	Recommended Dilution
<i>For use in automatic washers</i>	0.25% - 0.75% solution
<i>Descaling equipment</i>	0.75% - 4.0% solution

Packaging

Packaging	Product Code	Case Qty
Acid Tray Wash EVO 1 gallon	418086/1	4
Acid Tray Wash EVO 5 gallon	418086/5	1
Acid Tray Wash EVO 55 gallon	418086/55	1

Active:	Glycolic Acid
Use:	Acidic non-foaming cleaner for automated tray washers or general farm use
Activity:	General cleaner, hard water scale, soap scum and mineral deposit remover
Application:	Automated tray washers or pressure washing equipment

Benefits

- Phosphate-free formulation
- Ideal for use through chemical proportioners
- Effectively cleans trays, racks, plastic flats and boxes in automatic washers
- Removes mineral scale, detergent residues and other build-up from surfaces
- Formulated using biodegradable detergents

