



BioSentry® BioPhene™ Spray Disinfectant

Phenolic Disinfectant

Product

BioSentry® BioPhene™ Spray Disinfectant is an EPA-registered aerosol disinfectant that effectively controls the common viruses and bacteria that challenge the agriculture industry. In addition to controlling viruses, bacteria, and fungi, BioPhene Spray Disinfectant also eliminates odor and controls mold and mildew. BioPhene Spray Disinfectant can be used in vehicles, on footwear, or anywhere quick disinfection may be needed which makes it a great addition to your biosecurity portfolio.

Active: O-Phenylphenol and Ethanol

Use: Service vehicles, animal transport vehicles, footwear, animal care facilities, livestock facilities, office quarters, animal quarantine areas, animal laboratories, and animal life science laboratories

Activity: Virucidal, Tuberculocidal, Fungicidal, Staphylocidal, Pseudomonacidal

Application: Aerosolized Spray

Safety: Danger – corrosive, causes irreversible eye damage, combustible

Benefits

- Provides fast and effective disinfection
- Effective against Influenza A and Influenza B Viruses
- Kills 99.9% of harmful bacteria
- Eliminates odors at the source
- Dries quickly leaving behind no residue

Packaging

Item No.	Description
490910	BioSentry BioPhene Spray Disinfectant 15.5 oz

Usage Rates

Task	Usage
<i>General Disinfection</i>	Apply spray from 6-8 inches away until surface is wet, allow to dry

Active Ingredients

O-Phenylphenol	0.1%
Ethanol	63.2%
Other Ingredients	36.7%
Total	100.00%



800-621-8829 (USA/Canada) or 859-254-1221
 inform@neogen.com | animalsafety.neogen.com